**Profitest TN50M. W02**  
**Multilingual**  
**(Nasal, Nasen-Rachen)**

WANTAI SARS-CoV-2 Ag Schnelltest (Kolloidales Gold)  
von Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.


**AT Nr.: AT496/20**  
**PZN: 17942189**  
**HSC common list (RAT-ID: 1485) Category A**

**MedRhein**


**Profitest TN50M. W02**  
**Multilingual**  
**(Nasal, Nasen-Rachen)**

WANTAI SARS-CoV-2 Ag Schnelltest  
(Kolloidales Gold)

Testkassette  
X 50

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# HSC common list Vom 10.Juni.2022

EUROPEAN COMMISSION  
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY

Public health  
**Health Security**

## EU HEALTH PREPAREDNESS

### **EU Common list of COVID-19 rapid antigen tests**

*and a list of mutually recognised COVID-19 laboratory based antigenic assays*

*Agreed by the Health Security Committee*

*Last update: 10 June 2022*

# 1. Introduction

## 1.1 The EU common list of COVID-19 rapid antigen tests

This document builds on the Council Recommendation of 21 January 2021<sup>1</sup>, for which EU Member States unanimously agreed to set a common framework for the use and validation of rapid antigen tests and the mutual recognition of COVID-19 test results across the EU. The Council Recommendation called for a common framework that consists of COVID-19 rapid antigen tests that carry CE marking and that meet defined minimum performance requirements. Moreover, the devices should have been validated through an independent study carried out in at least one EU Member State, their use should be considered appropriate in the context of the COVID-19 pandemic, and their use should be in line with countries' national testing strategies.

On 17 February 2021, and on the basis of the Council Recommendation of 21 January 2021, the Health Security Committee agreed on a common list of COVID-19 rapid antigen tests (see Annex I). Since then, the EU common list has been regularly updated, taking into account the results of new validation studies and new rapid antigen test devices entering the market, as well as epidemiological developments and the emergence of SARS-CoV-2 variants.

The devices included in the EU common list of COVID-19 rapid antigen tests meet the criteria as defined by the Council Recommendation of 21 January 2021 as well as further criteria agreed by the Health Security Committee on 21 September 2021. A dedicated Technical Working Group on COVID-19 diagnostic tests<sup>2</sup> was set up by Health Security Committee with the objective to assess proposals submitted by countries and manufacturers for new devices to be included in the EU common list against these agreed criteria.

The EU common list of COVID-19 rapid antigen tests has been split up in two categories:

- Category A: COVID-19 rapid antigen tests evaluated by prospective clinical field studies; and
- Category B: COVID-19 rapid antigen tests evaluated by retrospective in vitro studies.

Since October 2021, the Health Security Committee has also agreed on a list of mutually recognised COVID-19 laboratory based antigenic assays (Annex II), which are meeting the same criteria as the COVID-19 rapid antigen tests.

## 1.2 EU Digital COVID Certificates

As determined by Regulation (EU) 2021/953<sup>3</sup>, the devices included in the EU common list of COVID-19 rapid antigen tests (Annex I) and carried out by health professionals or by skilled testing personnel, can be used by EU Member States to issue EU Digital COVID test

---

<sup>1</sup> Council Recommendation of 21 January 2021 on a common framework for the use and validation of rapid antigen tests and the mutual recognition of COVID-19 test results in the EU (OJ C 24, 22.1.2021, p.1).

<sup>2</sup> [https://ec.europa.eu/health/health-security-and-infectious-diseases/crisis-management/covid-19-diagnostic-tests\\_en](https://ec.europa.eu/health/health-security-and-infectious-diseases/crisis-management/covid-19-diagnostic-tests_en).

<sup>3</sup> Regulation (EU) 2021/953 of the European Parliament and of the Council of 14 June 2021 on a framework for the issuance, verification and acceptance of interoperable COVID-19 vaccination, test and recovery certificates (EU Digital COVID Certificate) to facilitate free movement during the COVID-19 pandemic (OJ L 211, 15.6.2021, p. 1–22).

# **strongly encouraged: rapid antigen tests included under Category A**

certificates. Moreover, as of 22 February 2022<sup>4</sup>, following a positive result of a COVID-19 rapid antigen test included in the EU common list and carried out by health professionals or by skilled testing personnel, it is possible for EU Member States to issue EU Digital COVID recovery certificates. These certificates may be issued retroactively, based on rapid antigen tests carried out from 1 October 2021. Both the EU Digital COVID test certificates and the EU Digital COVID recovery certificates should include specific data fields, as set out in the Annex to Regulation (EU) 2021/953.

As stipulated in point 15 of the Council Recommendation of 21 January 2021, Member States will agree on a selection of rapid antigen tests of which they will mutually recognise the test results for public health measures. The Health Security Committee agrees that, considering that *all* of the rapid antigen tests included in the EU common list are eligible for issuing EU Digital COVID test and recovery certificates, the entire list is considered to consist of rapid antigen tests of which Member States mutually recognise the test results for public health measures.

EU Member States are strongly encouraged to use, in particular, rapid antigen tests included under Category A of the EU common list for the issuance of EU Digital COVID certificates.

Secondly, EU Member States should pay particular attention to the issuance of EU Digital COVID recovery certificates based on the result of devices listed under Category B and that have solely been evaluated by the Paul-Ehrlich-Institut (PEI) in Germany, as only the sensitivity of these rapid antigen tests has been evaluated. Thirdly, EU Member States are strongly encouraged to ensure that only test results from the evaluated specimen type(s) are used for the issuance of EU Digital COVID test and recovery certificates.

It is important to note that results of COVID-19 laboratory based antigenic assays included in Annex II of this document are, at the moment, not eligible for the issuance of EU Digital COVID certificates.

## **2. The EU common list of COVID-19 rapid antigen tests**

### **2.1 Criteria to be met**

Based on a proposal by their Technical Working Group and taking into account the criteria presented by the Council Recommendation of 21 January 2021, the Health Security Committee agreed on 21 September 2021 on an updated scope and definitions of the EU common list of COVID-19 rapid antigen tests, as well as the criteria to be met by the applications for new rapid antigen tests to be possibly included in the EU common list.

The scope, definitions and criteria agreed on 21 September 2021 and as set out in this chapter have been applied to all proposals received by manufacturers and countries since the EU common list of COVID-19 rapid antigen tests was first agreed on 17 February 2021.

---

<sup>4</sup> Commission Delegated Regulation (EU) 2022/256 of 22 February 2022 amending Regulation (EU) 2021/953 of the European Parliament and of the Council as regards the issuance of certificates of recovery based on rapid antigen tests (OJ L 42, 23.2.2022, p. 4–8).

## ANNEX I: EU common list of COVID-19 rapid antigen tests <sup>12</sup>

**Disclaimer:** The Technical Working Group strongly recommends that rapid antigen tests are primarily used for preliminary testing for SARS-CoV-2 infection in symptomatic patients, and notes that rapid antigen tests should in particular be used in the specific contexts and circumstances referred to by the Commission Recommendation (EU) 2020/1743 and the updated technical report by ECDC on 26 October 2021. The content of the EU common list is based on the clinical performance data and information that is available at this moment in time. The Medical Device Coordination Group Guidance on performance evaluation of SARS-CoV-2 in vitro diagnostic medical devices<sup>13</sup>, envisaged to form the basis for common specifications to be adopted according to Article 9 of Regulation (EU) 2017/746, has been taken into consideration in this regard.

### Category A: COVID-19 rapid antigen tests evaluated by prospective clinical field studies

EU Member States are strongly encouraged to use, in particular, the rapid antigen tests included under “Category A” of the EU common list for the issuance of EU Digital COVID certificates. The clinical performance of these devices – for the specimen type as indicated in the corresponding column - has been evaluated by (at least) one prospective clinical field study meeting the criteria and definitions as agreed by the Health Security Committee on 21 September 2021.

EU Member States are strongly encouraged to ensure that only test results from the evaluated specimen type(s) are used to issue EU Digital COVID certificates.

Devices highlighted in blue are identical in design and construction but are, for example, branded or distributed under a different name. The results of validation studies may be transferred between devices that are identical in design and construction.

Device ID # <sup>14</sup>	Name of submitting company (and role) <sup>15</sup>	Commercial name of the device <sup>15</sup>	Clinical performance of the device <i>As evaluated by independent validation studies, meeting the agreed criteria</i>	Evaluated specimen type(s) <i>Eligible for issuing EU Digital COVID certificates</i>	Other specimen type(s) <sup>15</sup> <i>Not evaluated</i>	SARS-CoV-2 Target protein	Included in the EU common list since D/M/Y
1833	AAZ-LMB	COVID-VIRO®	<b>Prospective clinical field study</b> Prospective study carried out in the “Centre Hospitalier d’Orléans” on nasopharyngeal swabs, simultaneously tested by RT PCR: sensitivity <7 days after onset of symptoms: 94.7%, specificity: 100%.	Nasopharyngeal	Nasal	Nucleocapsid protein	10/05/2021
1232	Abbott Rapid Diagnostics (manufacturer)	Panbio™ COVID-19 Ag Rapid Test	<b>Prospective clinical field study</b> Study enrolling 1367 and 208 subjects in Utrecht (NL) and Aruba, respectively. NP swabs.	Nasopharyngeal	Nasal	Nucleocapsid protein	17/02/2021

<sup>12</sup> This is the list of rapid antigen tests as referred to in Article 3 of the Regulation (EU) 2021/953 of the European Parliament and of the Council of 14 June 2021 on a framework for the issuance, verification and acceptance of interoperable COVID-19 vaccination, test and recovery certificates (EU Digital COVID Certificate) to facilitate free movement during the COVID-19 pandemic, OJ L 211, 15.6.2021, p. 1–22.

<sup>13</sup> [https://ec.europa.eu/health/sites/default/files/md\\_sector/docs/mdcg\\_2021-21\\_en.pdf](https://ec.europa.eu/health/sites/default/files/md_sector/docs/mdcg_2021-21_en.pdf).

<sup>14</sup> As registered in and used by the JRC database; see: <https://covid-19-diagnostics.jrc.ec.europa.eu/>.

<sup>15</sup> Identical to what is included in the Instructions For Use (IFU) and/or labelling of the rapid antigen test.

Device ID # <sup>14</sup>	Name of submitting company (and role) <sup>15</sup>	Commercial name of the device <sup>15</sup>	Clinical performance of the device <i>As evaluated by independent validation studies, meeting the agreed criteria</i>	Evaluated specimen type(s) <i>Eligible for issuing EU Digital COVID certificates</i>	Other specimen type(s) <sup>15</sup> <i>Not evaluated</i>	SARS-CoV-2 Target protein	Included in the EU common list since <i>D/M/Y</i>
1778	Beijing Kewei Clinical Diagnostic Reagent Inc	COVID19 Antigen Rapid Test Kit	<p><b>Prospective clinical field study</b> Study at the General Hospital Jesenice in Slovenia: 103 RT-PCR positives and 450 RT-PCR negative subjects; symptomatic patients only. Overall sensitivity: 91.26%; specificity: 99.33%.</p> <p><b>Retrospective in vitro study</b> Positive evaluation by Paul-Ehrlich-Institut (PEI) in Germany: Sensitivity of 100% at Ct ≤ 25; <i>Manufacturer specificity of 100%.</i></p>	Nasal	-	Nucleocapsid protein	21/12/2021
1485	Beijing Wantai Biological Pharmacy Enterprise Co., Ltd	Wantai SARS-CoV-2 Ag Rapid Test (colloidal gold)	<p><b>Prospective clinical field study</b> Independent prospective study by Public Health Institute Ostrava (Czechia), including nasopharyngeal swabs from unselected symptomatic and asymptomatic participants. Sensitivity 80.6%, specificity 98.5% on 155 positive and 325 negative samples against RT-PCR (N total = 480).</p> <p><b>Retrospective in vitro study</b> Positive evaluation by Paul-Ehrlich-Institut (PEI) in Germany: Sensitivity of 100% at Ct ≤ 25; <i>Manufacturer specificity of 98.2%.</i></p>	Nasopharyngeal	Nasal	Nucleocapsid protein	14/07/2021
2519	BIOLAN HEALTH, S.L.	COVID-19 Antigen Rapid Test (Colloidal Gold Method)	<p><b>Prospective clinical field study</b> Prospective study performed in Hospital Universitario de Cruces (Spain). Nasal specimen, 314 negative samples and 116 positive samples. Sensitivity 98.1% at Ct&lt;25; overall sensitivity 81%; specificity 98.1%.</p>	Nasal	-	Nucleocapsid protein	04/03/2022
2035	BioMaxima SA	SARS-CoV-2 Ag Rapid Test	<p><b>Prospective clinical field study</b> Study in Poland performed on 480 samples of NP swabs taken from symptomatic patients and from asymptomatic people in contact with an infected person. Positive results were obtained in 205 patients and in the molecular test 213 people. Negative results were obtained in 275 people and in the molecular test 267 people. Diagnostic sensitivity: 93.43% (95% CI: 91.61%~97.19%) and diagnostic specificity: 97.75% (95% CI: 93.74%~98.92%). <a href="#">Source</a>.</p> <p><b>Retrospective in vitro study</b></p>	Nasopharyngeal	-	Nucleocapsid protein	23/07/2021

Date: June 8<sup>th</sup>, 2022

## Statement

To whom it may concern,

**We, Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.** registered at **No. 31 Kexueyuan Road, Changping District, Beijing 102206, China**, as one of the leading Chinese enterprises specialized great progress in both product development and market presence. We've been committed to providing professional IVD solution to customers over the world.

All viruses, including SARS-CoV-2, the virus that cause COVID-19, change over time. Some changes may effect the virus's properties, such as how easily it spreads, the associated disease severity, or the performance of vaccines, therapeutic medicines, diagnostic tools, or other public health and social measures. The identified SARS-CoV-2 variants are as follows:

WHO label	Pango lineage	GISAID clade/lineage	Nextstrain clade
Alpha	B.1.1.7	GRY	20I(V1)
Beta	B.1.351	GH/501Y.V2	20H(V2)
Gamma	P.1	GR/501Y.V3	20J(V3)
Delta	B.1.617.2	G/478K.V1	21A
Lambda	C.37	GR/452Q.V1	21G
Eta	B.1.525	G/484K.V3	21D
Iota	B.1.526	GH/253G.V1	21F
Kappa	B.1.617.1	G/452R.V3	21B
Mu	B.1.621	GH	21H
Omicron*	B.1.1.529	GR/484A	21K, 21L, 21M, 22A, 22B, 22C

\*Includes BA.1, BA.2, BA.3, BA.4, BA.5 and descendent lineages.

hereby we declare that, the above-listed variants can be tested by WANTAI SARS-CoV-2 Ag Rapid Test (Colloidal Gold), code WJ-2950, and WANTAI SARS-CoV-2 RT-PCR Rapid test, code WS-2248.

Sincerely,

Beijing Wantai Biological Pharmacy Enterprise Co., Ltd.

# Declaration of Conformity

## CE

### Manufacturer:

BEIJING WANTAI BIOLOGICAL PHARMACY ENTERPRISE CO.,LTD, located at No. 31  
Kexueyuan Road, Changping District, Beijing, 102206, P.R. of China  
TEL +86 10 59528888 FAX +86 10 89705849

### European Representative:

Qarad BV Pas 257, 2440 Geel, Belgium

### Products Name:

Wantai SARS-CoV-2 Ag Rapid Test (Colloidal Gold)

Product codes: WJ-2950 ,

Registration number: BE-CA01 /1-90054-00008-IVD

### Classification Under IVDD:

Not listed in Annex II

### Conformity assessment route: Annex III

We hereby declare that the product mentioned above meets the provisions of the European Directive 98/79/EC for in vitro Diagnostic Medical Devices. All supporting documentations are retained at the premises of the manufacturer.

### General applicable directive:

Directive 98/79/EC of European Parliament and of the Council of 27 October 1998 on in vitro Diagnostic Medical Devices.

### Standards applied:

EN ISO 18113-1:2011'EN ISO18113-2:2011 EN ISO 15223-1:2016 EN 13612: 2002 EN ISO 23640:2015'EN 13641 :2002'EN ISO 14971 :2012'EN ISO 13485:2016.

Ms. Zhao Lingzhi (Vice General Manager)

Beijing, October 12, 2021