



ABOUT US

Xiamen Wiz Biotech Co., Ltd. is a high-tech biomedical company, devoted to the field of rapid diagnostic reagents and instrument. Wiz is located in Xiamen, China, established in July 2013, which is a listed enterprises in National Equities and Quotations. As a technology-driven company that prides itself on its scientific excellence, Wiz focused on technological innovation and roduct innovation with 15 authorized patents, each of our products embodies the creativity of our excellent scientists, who are working hard to continuously bring novel products to the China and world markets. SARS-CoV-2 rapid detection series of products, including SARS-CoV-2 antigen, antibody and influenza differential detection, are produced by Wizbiotech, which can easily and quickly carry out screening of SARS-CoV-2 in a large number of people, adding strength to epidemic prevention and control.

Xiamen Wiz Biotech Co., Ltd.

COVID-19 Antigen Self-Testing Kit



Self-testing




SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)

- ✓ Anterior nasal swab specimen, noninvasive
- ✓ Very simple to use
- ✓ Convenient, no instrument required
- ✓ Rapid, results within 15~20 minutes
- ✓ Cost-efficient

<http://en.wizbiotech.com/>

CE 1434

QUICK REFERENCE GUIDE



CLINICAL PERFORMANCE

Test Results	Reference PCR Results		
	Positive	Negative	Total
Positive	113	0	113
Negative	2	456	458
Total	115	456	571

Sensitivity: **98.26%** (95% C.I. 93.86%~99.79%)
 Positive Predictive Value: **100%** (95% C.I. 96.79%~100.00%)
 Overall Percent Agreement: **99.65%** (95% C.I. 98.74~99.96%)

Specificity: **100.00%** (95% C.I. 99.19%~100.00%)
 Negativity Predictive Value: **99.56%** (95% C.I. 98.43%~99.95%)

PACKING INFORMATION

Specification	Box size(mm)	Carton size(mm)	Carton G.W.(kg)	Boxes/Carton	Tests/Carton
1Test	150*70*15mm	615*315*370mm	14kg	400	400
2Tests	150*65*30mm	675*320*320mm	11kg	200	400
5Tests	170*80*40mm	420*360*340mm	9kg	80	400
10Tests	150*100*80mm	420*320*420mm	7kg	40	400
20Tests	210*140*70mm	440*300*370mm	7kg	20	400

CERTIFICATE



CERTIFICATE

EC Certificate No. 1434-IVDD-492/2021
EC Design-examination
Directive 98/79/EC concerning
in vitro diagnostic medical devices

Polish Centre for Testing and Certification certifies
that manufactured by:

Xiamen Wiz Biotech Co., Ltd.,
3-4 Floor, NO.16 Building, Bio-medical Workshop,
2030 Wengjiao Xi Road, Haicang District, Xiamen City,
Fujian Province, 361026, P.R. China

in vitro diagnostic medical devices
for self-testing

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)

51332801, 51332802, 51332803, 51332804

in terms of design documentation, comply with requirements
of Annex III (Section 6) to Directive 98/79/EC (as amended)
implemented into Polish law,
as evidenced by the audit conducted by the PCBC

Validity of the Certificate: from 22.11.2021 to 27.05.2024

The date of issue of the Certificate: 22.11.2021

The date of the first issue of the Certificate: 22.11.2021



Issued under the Contract No. MD-77/2021
Application No: 130/2021
Certificate bears the qualified signature.
Warsaw, 22/11/2021
Module A1

Vice-President

POUSH CENTRE FOR TESTING AND CERTIFICATION 02-844 Warsaw, 469 Puławska Street, tel. +48 22 46 45 200, e-mail: pcbc@pcbc.gov.pl

Certificate CN1942110

The management system of
Xiamen Wiz Biotech Co., Ltd.
3-4 Floor, NO. 16 Building, Bio-medical Workshop, 2030 Wengjiao Xi Road,
Haicang District, Xiamen City, Fujian Province, 361026, P.R. China

has been assessed and certified as meeting the requirements of

ISO 13485:2016
EN ISO 13485:2016

For the following activities

The scope of registration appears on page 2 of this certificate.

This certificate is valid from 26 August 2020 until 3 September 2022
and remains valid subject to satisfactory surveillance audits.
Re certification audit due before 26 July 2022
Issue 2, Certified since 4 September 2019

This is a multi-site certification.
Additional site details are listed on the subsequent page.

Authorized by

SGS United Kingdom Ltd
Rosemead Business Park, Elmestree Park, Chesham, CH85 5EN, UK
T +44 (0)191 250-6000 F +44 (0)191 250-6000 www.sgs.com

HC SGS 13485 2016 0118 M2
Page 1 of 2

This document is issued by the Company subject to the General Conditions of Certification Services available at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the importance of clarity, interpretation and confidentiality of these conditions. The validity of this document may be verified at www.sgs.com/verification or by contacting your local SGS office. Any modification, alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

Certificate CN1942110, continued

Xiamen Wiz Biotech Co., Ltd.

ISO 13485:2016
EN ISO 13485:2016

Issue 2

Detailed scope

Design and manufacture of in vitro diagnostic medical devices including hormone marker test kit, kidney disease marker test kit, tumor marker test kit, infectious disease marker test kit, gastrointestinal inflammation marker test kit, cardiac marker test kit, immunochromatography analyzer.

Additional facilities

4 Floor, NO. B12 Building, Bio-medical Workshop, 2072 Wengjiao Xi Road, Haicang District, Xiamen City, Fujian Province, 361026, P.R. China

Authorized by

UKAS MANAGEMENT SYSTEMS
0005

SGS United Kingdom Ltd
Rosemead Business Park, Elmestree Park, Chesham, CH85 5EN, UK
T +44 (0)191 250-6000 F +44 (0)191 250-6000 www.sgs.com

HC SGS 13485 2016 0118 M2
Page 2 of 2

This document is issued by the Company subject to the General Conditions of Certification Services available at www.sgs.com/terms_and_conditions.htm. Attention is drawn to the importance of clarity, interpretation and confidentiality of these conditions. The validity of this document may be verified at www.sgs.com/verification or by contacting your local SGS office. Any modification, alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

FACTORY

Our one hundred thousand grade clean and dust-free workshop is operated strictly under ISO 13485:2016 and GMP guidelines. With the automatic production line, efficient production process, strict quality control, we always produce the products timely to meet customers' request.



EN

EN

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold)

Instructions for Use

For use at home self-test or Non-professional
For use with nasal cavity (anterior nasal) swab specimen
For In Vitro Diagnostic Use Only

INTENDED USE

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) is intended for the qualitative detection of SARS-CoV-2 Antigen (Nucleocapsid protein) which is in nasal cavity (anterior nasal) swab specimen from individuals with suspected COVID-19 infection. The test kit is intended for self-test or home test.

SUMMARY

COVID-19 is a acute respiratory infectious disease. People are prone to infection generally. Currently, the main cause of COVID-19 infection is contact with someone who is already infected with the SARS-CoV-2, and asymptomatic infected people may also transmit the virus. Studies have shown that symptoms of infection generally appear within 1-6 days, with most occurring within 3 to 7 days after infection. The main symptoms are fever, fatigue, loss of smell or taste, and a dry cough. In some cases, a stuffy nose, runny nose, muscle pain and diarrhea can also occur.

PRINCIPLE THE DETECTION

SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) adopts immune lateral chromatography technology. When the test specimen contains SARS-CoV-2, the SARS-CoV-2 Antigen will react with the antibody coated on the test line (T) to appear a red band in the test line (T) area, when the content of the SARS-CoV-2 in the test specimen is too low or does not exist, the test line (T) area does not appear a red band. Regardless of whether the test specimen contains the SARS-CoV-2, a red band will appear in the area of the quality control line (C), which is the basis for judging whether the test is effective.

MAIN KIT COMPONENTS

Table with 4 columns: Catalogue number, 51332801, 51332802, 51332803, 51332804. Rows include Test kit, Test kit, Disposable swab, and Biowaste Waste Bag.

Materials required but not provided: Timer

PRECAUTIONS

- 1. For In Vitro Diagnostic Use Only.
2. For use with nasal cavity (anterior nasal) swab specimen.
3. For the detection of presence from SARS-CoV-2 only, not for other viruses or pathogens.
4. The person under 18 years of age should be tested with the assistance of the legal guardian or authorized person.
5. Keep the test kit or kit components out of the reach of children and pets before and after use.
6. Test card packaging contains desiccant, it is forbidden to eat.
7. The sample extraction Solution in extraction tube contains chemical components. Direct contact should be avoided and eat prohibited. If the solution contacts the skin or mucosa or eye, please flush with copious amount of water. Please contact your family doctor or professional or seek medical advice if necessary.
8. Use of gloves and other protective equipment are recommended when conducting testing.
9. The test kits should be stored according to storage conditions required in Instructions for Use. It is forbidden to use the test kit that is not stored as required.
10. Don't use the test kit beyond its expiration date.
11. Don't use kit components that have been opened or changed.
12. Leave test card sealed in its foil pouch until just before use. Don't use if pouch is damaged or open.
13. Disposable swabs are sterile products. Don't use if the swab packaging is damaged or opened.
14. The use of the swab should strictly follow instruction for use; otherwise it may cause nasal cavity bleeding, swab rupture and retraction and other risks.
15. Don't immerse the swab in the extraction solution or other liquid before collect sample with Disposable swab.
16. Don't touch swab soft tip when handling the swab sample.
17. Proper sample collection and handling are crucial for correct results.
18. Don't mix components from different kits.

- 19. All kit components are single use items. Don't use with multiple specimens. Don't reuse the used test kit or kit components.
20. Before deciding to implement relevant treatment or management decisions, it is recommended to conduct communication with family doctors or professionals. Don't take medicine privately or any action that is dangerous yourself or others.
21. The used test kit components and samples can be put into plastic bags together with ordinary household waste. If the test result is positive, you should carefully dispose the relevant waste components and samples, and thoroughly clean and disinfect the working surface to ensure hygiene. If there are special requires on waste in local laws and regulations, you must strictly observe it.
22. In view of the global epidemic of COVID-19, all actions should comply with the current measures and regulations of your country/region, to implement prevention and control measures scientifically, and protect yourself and others effectively.
23. Do not eat, drink or smoke in the area where handling specimens or test kits.

STORAGE CONDITIONS AND SHELF LIFE

The test kit should be stored conditions of 2°C~30°C, dry and out of direct sunlight (Don't freeze the kit or its components).
The shelf life of the kit is 12 months.
The test card should be used within 60 minutes after opening the aluminum foil bag.
For the kit expiration date, please refer to the product label.

FREQUENTLY ASKED QUESTIONS (FAQ)

- 1. What are the known or potential benefits of product testing?
2. What are the known or potential risks in product testing?
3. What factors will affect the test results? What should I pay attention to?
4. What's the difference between Antigen and Molecular test?
5. I have taken the test, but I don't see the control line (C). What should I do?
6. I have taken the test, but I don't see the test line (T). What should I do?
7. I have taken the test, but I don't see the test line (T) and control line (C). What should I do?
8. I have taken the test, but I don't see the test line (T) and control line (C). What should I do?
9. If the test result is positive, what should I do?
10. If the test result is negative, what should I do?
11. If the test result is positive, what should I do?
12. If the test result is negative, what should I do?
13. If the test result is positive, what should I do?
14. If the test result is negative, what should I do?
15. If the test result is positive, what should I do?
16. If the test result is negative, what should I do?
17. If the test result is positive, what should I do?
18. If the test result is negative, what should I do?
19. If the test result is positive, what should I do?
20. If the test result is negative, what should I do?

POSITIVE VALUE

Positive value/limit of detection: 1.7x10^7 TCID50/mL.
Select the confirmed activated SARS-CoV-2 medium (concentration 1.8x10^7 TCID50/mL) and use gradient dilution method to find out the virus medium which is reach the critical value of the detection. That is repeating the test for 20 times and the test result is positive for at least 19 times.

CLINICAL PERFORMANCE

The clinical performance characteristic of the SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) was established with 571 anterior nasal swab specimens. Specimens were collected and tested according to the requirements of the Instructions for Use. The storage, transportation and detection of samples after collection met the relevant requirements of the Instructions for Use. At the same time, the nasopharyngeal swab specimens of the same donor were detected by nucleic acid detection reagent (Virusay ePCR SARS-CoV-2). WIZ's SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) was compared with reference reagent PCR results.

Clinical performance of the WIZ'S SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) compared to RT-PCR

Table with 3 columns: WIZ Results, Reference PCR Results, Total. Rows include Positive, Negative, Total.

Sensitivity: 98.26% (95%CI: 93.86%/99.79%)
Specificity: 100.00% (95%CI: 99.19%/100.00%)
Positive Predictive Value: 100% (95%CI: 96.79%/100.00%)
Negative Predictive Value: 99.56% (95%CI: 98.43%/99.95%)
Overall Percent Agreement: 99.63% (95%CI: 98.74%/99.96%)

PERFORMANCE CHARACTERISTICS

- Using antisera reference for testing, the results meet the requirements of enterprise reference.
Cross reaction.

Table with 3 columns: Name (Microorganism), Concentration, Test result. Lists various pathogens like Influenza B Virus, Influenza A H1N1, SARS Coronavirus, etc.

Table with 3 columns: Name, Concentration, Test result. Lists various pathogens like Bacillus pumilus, Clostridium botulinum, Escherichia coli, etc.

Interference Substances

Table with 4 columns: Interference Substances, Concentration, Negative interference results, Positive interference results. Lists substances like Milk, Urine, Blood, etc.

Hook effect: Within the concentration of 3.4x10^7 TCID50/mL, cell culture medium of SARS-CoV-2 Antigen, the test results of this product showed no Hook effect.
The product was established by using internal reference, there were no differences observed within run.

LITERATURE REFERENCES

- [1] MPA. The Technical Key Points for Coronavirus (COVID-19) Antigen-antibody Detection Reagent Registration (2020).
[2] Xu Chu, Li Fan. Analysis on the Risk Management of In Vitro Diagnostic Reagents. China Medical Device Information, 2020, 36(18): 19-20.
[3] Wu, Jheng, Meng, H. Immunocolloidal Gold Technology: Advances and Application(s). Chinese Agricultural Science Bulletin, 2019, 35(12): 146-151.
[4] Yoon, J. Rapid Test of Immunocolloidal Gold with Membrane as Solid Phase Carrier(s). Progress in Microbiology and Immunology, 2003, 31(1): 74-78.


TEST PROCEDURE

It is important to read the Instructions for Use carefully and follow the steps in the correct order. It takes about 15 minutes to prepare before each test, and the results can be obtained after 30 minutes.

Table with 4 columns: Test Procedure, Preparatory work, Sample collecting, Sample processing, Sample testing.

I Preparatory work

- 1. Please use the test kit at room temperature (15°C~30°C). If the test kit was previously stored in a cool place (temperature below 15°C), please balance it at 15°C~30°C for 30 minutes before use.
2. Prepare a Timer (such as watch, clock), tissue, hand sanitizer/soap and warm water.
3. Please read the Instruction for Use carefully.
4. Wash hands thoroughly (at least 20 seconds) with soap and warm water /hand sanitizer. This step ensures that the kit will not be contaminated, then dry your hands.



Warning!


If the sealed package is damaged, don't use it. Please replace it with a new kit. If you have cough symptoms, get tested in private.

II Sample collecting

- 1. Take out the sample extraction tube, unscrew the tube cap of the extraction tube.
2. Press the swab soft head tightly against the inner wall of the extraction tube, and rotate the swab clockwise or anticlockwise for about 10 times.
3. Squeeze the swab head along the inner wall of the sample extraction tube under the liquid in the tube as much as possible, take out the swab.
4. Tighten the tube cap for standby.

III Sample processing

- 1. Take out the extraction tube, dip the swab soft head into the extraction solution and immerse it in the liquid.
2. Press the swab soft head tightly against the inner wall of the extraction tube, and rotate the swab clockwise or anticlockwise for about 10 times.
3. Squeeze the swab head along the inner wall of the sample extraction tube under the liquid in the tube as much as possible, take out the swab.
4. Tighten the tube cap for standby.




Warning!

Do not insert too far into the nostril as this may lead to nasal cavity bleeding or swab rupture and other risks. Select the confirmed activated SARS-CoV-2 medium (concentration 1.8x10^7 TCID50/mL) and use gradient dilution method to find out the virus medium which is reach the critical value of the detection. That is repeating the test for 20 times and the test result is positive for at least 19 times.

IV Sample testing

- 1. Tear off the aluminum foil bag, take out the test card and place it horizontally on the test desk.
Warning! The platform should be in a horizontal and stable state, and tilt and shake are strictly prohibited.
2. Unplug the adding sample hole cover of the extraction tube.
3. Gently squeeze the extraction tube, and drop 2 drops liquid vertically into the sample well of the test card.
Warning! The existence of bubbles in the extraction tube may lead to the wrong sample volume and inaccurate test results. If there are bubbles in the extraction tube, gently shake the extraction tube to squeeze out part of the liquid so as to remove the bubbles.
4. Start timing, read the test results at 15 minutes. Do not read result before 15 minutes or after 30 minutes.




INTERPRETATION OF THE RESULTS

Negative: The quality control line (C line) appears red band, while the test line (T line) doesn't appear red band.
Positive: Both the quality control line (C line) and the test line (T line) appear red bands.

Caution!

Regardless of the shades of color within the specified detection time, the result should be judged as positive result.



Warning! Invalid results may indicate that you did not strictly follow the operation steps of instruction for use. Please read the instruction for use carefully again and select a new kit for testing again.

PROCESSING OF TEST RESULTS

- If the test result is positive: COVID-19 infection is present currently. Contact your doctor or local health department immediately. Comply with the self-quarantine requirement and protection guidelines in your area. Carry out PCR test for confirmation.
If the test result is negative: You need to continue to take measures related to contact with others and self-protection. Even if your test result is negative, it is also possible that there is an infection. If you still suspect, please repeat the test after 1-2 days. On account of the coronavirus cannot be accurately detected at every stages of infection.
If the test result is invalid: It may be caused by incorrect operation in the detection process. Please repeat the test. If the test result is still invalid, please contact your doctor or COVID-19 testing center.

QUALITY CONTROL

The quality control line is a key point of test kit and is used to control the procedure. The quality control line appears, which indicates that the test has been performed correctly and the test kit has reached SYMBOLS.

Table with 6 columns: Symbol, Used for, Symbol, Used for, Symbol, Used for. Lists symbols for Control substance, Validity date, Calibration, Test kit, Test device, and Calibration date.

Xiamen Wize Biotech Co., Ltd.
The quality control line (C line) appears red band, while the test line (T line) doesn't appear red band.
Positive: Both the quality control line (C line) and the test line (T line) appear red bands.

Caution! Regardless of the shades of color within the specified detection time, the result should be judged as positive result.

Xiamen Wize Biotech Co., Ltd.
The quality control line (C line) appears red band, while the test line (T line) doesn't appear red band.
Positive: Both the quality control line (C line) and the test line (T line) appear red bands.

Supplier of Disposable Sterile Swab
Zhejiang Gongdong Medical Technology Co., Ltd.
CE 0123 (According to Directive 93/42/EEC)

Address: No. 10 Beiyuan Ave., Huangyan, 318020 Taizhou, Zhejiang, PEOPLE'S REPUBLIC OF CHINA
CEP Name: Shanghai International Holding Corp GmbH(Europe)
CEP Address: Effelstrasse 80,20537 Hamburg,Germany

CE 1434
CE 03.001.2021.10.156p